

*for those little ones  
in need of some peace & quiet...*

# MiniMuffs®



## NEONATAL NOISE ATTENUATORS

**DESIGNED  
SPECIFICALLY FOR  
PREMATURE BABIES,  
NATUS MINIMUFFS  
NOISE ATTENUATORS  
PROTECT THEIR  
SENSITIVE EARS,  
AND PROVIDE A  
SAFE ENVIRONMENT  
FOR HEALTHY  
DEVELOPMENT.**

### Provides a “quiet hour” for premature babies in the NICU

- Helps meet American Academy of Pediatrics guidelines for safe noise levels in the NICU<sup>1</sup>
- Reduces sound levels by at least 7 decibels (dB)<sup>1</sup>
- Results in a reduction of sound pressure level of more than 50%<sup>1</sup>

### Ideal during special patient care procedures and transport

- Protects the sensitive ears of babies under ventilation
- Decreases the MRI noise levels when used in conjunction with earplugs<sup>1</sup>
- Reduces disruptive noise associated with transfer of infants to different hospital units, and provides quieter environment during transport in emergency vehicles

### Safe and Comfortable

- Fits comfortably around the baby’s ears with a soft foam oval-shaped design, and gentle hydrogel adhesive
- Meets infection control standards through single-patient use design

<sup>1</sup>References on file at Natus

## The MiniMuffs Neonatal Noise Attenuators

### Pack includes:

- 36 individually wrapped pairs
- Complete Instructions for Use



### Specifications:

#### Nominal Dimensions:

- Inner width: 1.125 in. (2.856 cm)
- Inner length: 1.375 in. (3.493 cm)
- Inner height: 0.31 in. (0.787 cm)

Reduced noise levels are clinically demonstrated to positively impact infant state<sup>1</sup>. By reducing noise levels, MiniMuffs noise attenuators result in potentially longer periods of sleep time, optimal oxygenation, and improved physiological stability.

**MiniMuffs**<sup>®</sup>  
N E O N A T A L N O I S E A T T E N U A T O R S

*giving newborns  
the quiet time they need.*

**natus**<sup>®</sup>  
SOLUTIONS FOR NEWBORN CARE

Natus Medical Inc. • San Carlos, CA USA  
Natus U.K. • Oxford, England  
www.natus.com  
customer\_service@natus.com  
+1 650-802-0400, 800-255-3901  
Fax: +1 650-802-0401



P/N 050312F