

Richmond Ultra Low Bed



The Richmond Ultra Low Bed has been designed specifically to provide the ultimate solution for nursing patients who may be at risk of falling out of bed. With an exceptional lowest height position of only 16cm, it ensures a safe sleeping surface for the patient, helping to reduce the risk of fall related injuries. To facilitate the varying requirements of the hospital environment, the Richmond Ultra Low bed comes fully equipped with enhanced features such as electric CPR, quick release head and foot boards, integral side rails and a central braking system as standard.



Technical Specifications

Product Code	E5-8500	Knee Break Angle	33°
Weight Limit	28.3 stone / 180 kgs	Castors	7.5cm
Overall Length	208.5cm	Actuators	LINAK
Overall Width: Without Side Rails With Side Rails	90cm 101cm	Back Rest	4000N
		Knee Break	4000N
		Raise Up	6000N
Mattress Platform Dimensions	198cm x 88cm	Total Bed Weight	120 kgs
Platform Height: Low Position High Position	16cm 70cm	Transportation Frame: Weight Dimensions	8 kgs 112cm x 41cm x 25cm
Side Rail Height	41cm	Rating	IP66
Side Rail Length	169cm	Warranty	5 years
Back Rest Angle	70°		

*All equipment should be used in conjunction with clinical judgement and assessments dependent on individual patient morphology. Essential nursing care is pivotal in providing the optimum pressure ulcer prevention plan.



FEATURES	BENEFITS
Ultra low height position of 16cm	Helps to reduce the risk of injury by providing a safe sleeping surface for patients who may be at risk of falling out of bed
Simple and easy to use patient handset	Increases patient independence by allowing the patient to adjust the position of the bed
Nurse lock out facility	Locks the electrical functions of the bed and prevents unauthorised setting changes
Auto contour function	Reduces shear and friction forces by allowing the back rest and knee break sections to adjust simultaneously
Electrically adjustable back rest and knee break	Bed can be easily adjusted to suit all comfort requirements of the patient, also aids moving and handling procedures
Electrically variable height	Allows the height of the bed to be adjusted to suit all nursing procedures
Manual leg raiser	Offers a range of different positions to aid relief and comfort
Removable ABS mattress platforms	Each platform can be individually removed to allow easy cleaning and aid disinfection protocol
Quick release head and foot boards	Can be removed easily to ensure quick access for emergency procedures, also allows for easy cleaning
Electric CPR function under single button operation	Ensures quick levelling of the bed to aid emergency procedures
Manual CPR levers in back rest section	Easy access to ensure quick levelling of the bed in an emergency
Integral quick release side rails	Provides extra security for the patient and ensures conformance to MHRA guidelines. Also allows quick access for emergency procedures
Battery back up as standard	Allows the bed to be lowered in the event of a power failure to assist with health & safety issues
Central braking system	Provides a secure braking system and ensures excellent steering and manoeuvrability of the bed
Crash mat available	Maximises patient safety by helping to reduce the risk of injury
Drip mounting facility	Supports all nursing procedures
Lifting pole facility	Aids patient transfers and facilitates moving and handling procedures
Transporter kit as standard	Easily manoeuvrable and compact for handling and storage
IP66 Ingress protection on all electrical components	Aids decontamination protocol and infection control
Compatible with all Park House mattresses	Provides a full pressure care package



Park House
Healthcare

Tel: 0845 0600 333

Fax: 0845 0600 334

Email: sales@parkhouse-hc.com

Web: www.parkhouse-hc.com