



IntuBrite Laryngoscopes

Laryngoscope Handles by IntuBrite

E-Max Handle

IntuBrite has completely redesigned the laryngoscope handle in order to deliver a better lift angle, a more positive gripping surface, a more reliable power cell, and a precision driving circuit for the dual-LED blade set-up. Constructed of solid aluminum and offering a robust, click-in, steel connecting head, the new handle has been receiving rave reviews for comfort and performance from early users.

Each IntuBrite EMS kit will be supplied with one handle; additional handles will be available for purchase individually for back-up or so that individuals may carry their own in their personal trauma bags.



E-Flex Handle



The E-flex handle was ergonomically designed with pediatric patients in mind. The e-flex offers a slightly shorter handle with deeper finger grips. The gentle arc in the handle provides a perfect angle for proper technique and mandible lift. The top thumb groove allows for a lighter grip and greater control while attending to small infants and children. Although originally designed as a pediatric handle, we have found that many practitioners in our trial studies like the control they feel and leverage attained with the e-flex handle while attending to children and adults.

E-Classic Handle

During the process of development and field trials, a few anesthesiologists and seasoned paramedics alike requested the option of a "straight" handle; The E-Classic is IntuBrite's version. The E-Classic handle retains our signature "pistol" finger grips and adds a subtler angle for an ergonomic mandible lift in a more familiar design. Like all of IntuBrite's innovative handles the e-classic is CNC Machine tooled from a solid piece of block aluminum and the blade clicks onto a stainless steel bar mount providing a stable and secure fit for the blade.

All IntuBrite handles come with a lifetime guarantee.



STAT-TIAKENI MEDICAL (Pty) Ltd.

Suite 247, Private Bag X7, Northriding, 2162 A.P.D. Industrial Park, Unit 41 & 42

Elsecar & Kelvin Streets, Kya Sand

Tel: (011) 708-6363 Fax: (011) 708-6364

e-mail: sales@statmedical.co.za website: www.statmedical.co.za

Laryngoscope Blades by IntuBrite

Machintosh and Miller Laryngoscope Blades

IntuBrite™ manufactures its uniquely designed laryngoscope blades from surgical-grade, stainless steel. These are custom-fitted with a pair of proprietary, UV and White LED bulbs (blacklight and whitelight), with wavelengths specifically “tuned” to deliver precisely the correct illumination which provides their unique anatomy-differentiating and almost glareless view.

Specially molded medical-grade plastic has been bonded to the surgical stainless material to form the “snap-in” base which locks each blade to the laryngoscope handle prior to use and provides power transition from the handle to the lights in the blade. IntuBrite's blades have been tested to 55lbs. of pressure which exceeds regulations and is far superior to any other disposable on the market today.

IntuBrite manufactures a selection of 10 of the most popular sizes and styles of laryngoscope blades:

- MacIntosh 1
- MacIntosh 2
- MacIntosh 3
- MacIntosh 4
- Miller 00
- Miller 0
- Miller 1
- Miller 2
- Miller 3
- Miller 4



Various combinations of these types and sizes may be provided to EMS and ER/OR professionals depending upon their departmental and/or personal clinical preferences.

STAT-TIAKENI MEDICAL (Pty) Ltd.

Suite 247, Private Bag X7, Northriding, 2162 A.P.D. Industrial Park, Unit 41 & 42
Elsecar & Kelvin Streets, Kya Sand

Tel: (011) 708-6363 Fax: (011) 708-6364

e-mail: sales@statmedical.co.za website: www.statmedical.co.za