




flexicare[®]
TOTAL QUALITY - TOTAL CARE

suction
product catalogue



SUCTION



SUCTION YANKAUERS

130 SUCTION YANKAUERS

SUCTION CATHETERS



SUCTION CATHETERS **131**

SUCTION TUBING

131 SUCTION TUBING

Suction Yankauers

Effective suction without aspiration of surrounding tissue, and a choice of with or without fingertip vacuum control. Transparent, with a smooth inner lumen for uninterrupted suctioning and rapid passage of material.



Rosebud Tip



Straight Tip

Continuous Suction

The lateral eyes positioned near the tip for continuous and easy suctioning.

Optional Vacuum Fingertip Control

Port at the distal end of the handle gives fingertip / thumb control over the vacuum level.

Occlusion Discernibility

Blockages or obstruction are easily seen through the transparent tube. The smooth inner lumen helps maintain the level of suction.

STERILE WITH VACUUM CONTROL	STERILE STANDARD		QTY
050-35-001	050-35-002	Sterile Suction Yankauer with Straight Tip	25
050-35-003	050-35-004	Sterile Suction Yankauer with Rosebud Tip	25



Suction Catheters

A smooth outer surface, rounded tip to reduce risk of injury and soft edges to the double sided eyelets positioned to prevent blockages, makes the sterile Suction Catheter ideal for Oral, Tracheal or Endobronchial suctioning. Available as Standard or with Vacuum Control.

Smooth Clear Main Tube

Smooth outer surface for easy insertion through all types of endotracheal and tracheostomy tubes.

Vacuum Control

Fingertip port gives the user control over the vacuum level.

STERILE STANDARD	STERILE WITH VACUUM CONTROL		LENGTH	Fr	OUTER DIAMETER MM	COLOUR	QTY
050-10-06	050-20-06	Sterile Suction Catheter	500mm	6	2.00	●	100
050-10-08	050-20-08	Sterile Suction Catheter	500mm	8	2.66	●	100
050-10-10	050-20-10	Sterile Suction Catheter	500mm	10	3.33	●	100
050-10-12	050-20-12	Sterile Suction Catheter	500mm	12	4.00	○	100
050-10-14	050-20-14	Sterile Suction Catheter	500mm	14	4.66	●	100
050-10-16	050-20-16	Sterile Suction Catheter	500mm	16	5.33	●	100
050-10-18	050-20-18	Sterile Suction Catheter	500mm	18	6.00	●	100




Suction Tubing

Strong wall construction makes the suction tubing collapse resistant under high negative pressure and is ideal for the removal of fluids during medical procedures and treatments.

Spare Male-Male Connector

Supplied sterile for connection to Female connectors, or to connect two tubes if extended length is required.

7MM ID ●	6MM ID ●	5MM ID ●		QTY
1020-25	1020-6		Sterile Suction Tube 2m	25
1024-25			Sterile Suction Tube 2.4m	25
1027-25			Sterile Suction Tube 2.7m	25
1030-25	1030-6	1030-5	Sterile Suction Tube 3m	25
1035-25			Sterile Suction Tube 3.5m	25
1037-25			Sterile Suction Tube 3.7m	25
10830-7		10830-5	Suction Tube 30m coil with bubble non conductive*	1



Discover more online. Visit our website for the very latest product releases, instructional videos, company news, event information and downloadable resources.

Subscribe to our social media channels to join the conversation and share your stories.

www.flexicare.com



Flexicare Group Limited

Cynon Valley Business Park • Mountain Ash • CF45 4ER • UK
+44 (0)1443 474647 • enquiries@flexicare.com

www.flexicare.com

FC7087A | 112022