



2021/2023

RT₂-2021 - RT₄-2023

BROCHURE



SCAN ME



TRACT



ZOLL M2 - BASIC WITH PACING

TENDER NUMBER	RT2-05-04
DESCRIPTION	<ul style="list-style-type: none">- DEFIBRILLATOR- MONITOR- 3/5 LEAD ECG- BUILT IN AED- LITHIUM RECHARGEABLE BATTERY- ZOLL RECTILINEAR BIPHASIC TECHNOLOGY- EXTERNAL PADDLES- CARRY CASE- CPR FEEDBACK- PACING



RT2-2021 PRICE LIST

Replacement Accessories

RT2-05-04 - Zoll M2 Defibrillator with Pacing and ECG

RT2-2021 ITEM NUMBER	DESCRIPTION	MASTER MATERIAL NUMBER
RT2-05-04	ZOLL M2 DEFIBRILLATOR WITH PACING	220484
RT2-05-04a	Reusable external adult/pediatric defibrillator paddles	222481
RT2-05-04b	3-Lead ECG patient cable	222582
RT2-05-04c	ECG electrodes (x400)	222483
RT2-05-04d	AED pads	222485
RT2-05-04f	Pacing electrodes	222485
RT2-05-04h	Electrode gel (x6: 250ml)	222486
RT2-05-04i	Printer paper	222487
RT2-05-04j	Dedicated trolley	222488
RT2-05-04kj	Internal paddles (Pair)	TBA



RT2-2021 CONTRACT



ZOLL M2 with SPO₂

TENDER NUMBER	RT2-05-05
DESCRIPTION	<ul style="list-style-type: none">- DEFIBRILLATOR- MONITOR- 3/5 LEAD ECG- BUILT IN AED- CARRY CASE- PACING- SPO₂- LITHIUM RECHARGEABLE BATTERY- EXTERNAL PADDLES- CPR FEEDBACK



RT2-2021 PRICELIST

Replacement Accessories RT2-05-05 - Zoll M2 Defibrillator with Pacing, ECG and SP02

RT2-2021 ITEM NUMBER	DESCRIPTION	MASTER MATERIAL NUMBER
RT2-05-05	ZOLL M2 DEFIBRILLATOR WITH PACING & SPO2	220485
RT2-05-05a	3-LEAD PATIENT CABLE	222382
RT2-05-05b	ECG ELECTRODES	222483
RT2-05-05c	AED pads	222484
RT2-05-05e	Electrode gel	222486
RT2-05-05f	Printer paper	222487
RT2-05-05g	SpO2 transducer	222491
RT2-05-05h	SpO2 cable	151988
RT2-05-05i	Dedicated trolley	222488
RT2-05-05j	Internal paddles	TBA



RT2-2021 CONTRACT



ZOLL R- SERIES WITH NIBP & SPO₂

TENDER NUMBER
DESCRIPTION

RT2-05-06

- DEFIBRILLATOR
- MONITOR
- REAL-CPR-HELP
- 3 LEAD ECG CABLE
- CODE READINESS TESTING SYSTEM
- ZOLL RECTILINEAR BIPHASIC TECHNOLOGY
- BUILT-IN TEST PORT
- INTEGRATED AC POWER
- PACING
- MASIMO PULSE OXIMETRY (SPO₂)
- NON-INVASIVE BLOOD PRESSURE (NIBP)
- LITHIUM-ION BATTERY PACK
- EXTERNAL PADDLES




RT2-2021 PRICE LIST

Replacement Accessories

RT2-05-06 - Zoll R-series Defibrillator with ECG SPO₂ and NIBP

RT2-2021 ITEM NUMBER	DESCRIPTION	MASTER MATERIAL NUMBER
RT2-05-06	ZOLL R-SERIES WITH NIBP & SPO ₂	220486
RT2-05-06a	3-LEAD PATIENT CABLE	222382
RT2-05-06b	ECG ELECTRODES (x 400)	222482
RT2-05-06c	AED pads	222483
RT2-05-06d	Cable for AED pads	222484
RT2-05-06e	Electrode gel (x6: 250ml)	222486
RT2-05-06f	Printer paper (x10 packs)	222487
RT2-05-06g	SpO ₂ sensor	152080
RT2-05-06h	SpO ₂ cable	151988
RT2-05-06i	BP hose	151713
RT2-05-06j	BP cuff (Adult/ Peads)	151713
RT2-05-06ak	Dedicated trolley	222488
RT2-05-06l	Internal paddles	TBA
RT2-05-06m	Reuseable external adult/pediatric Defibrillator paddles	222481

RT4-2023 CONTRACT



SAM IO NEEDLES

TENDER NUMBERS: RT4-06-044 Adult
 RT4-06-045 Paediatric

DESCRIPTION:



These needles are a stylet/catheter combination **meant to be used with the SAM IO DRIVER.**





The snap-fit connection keeps the needle in place while the manually actuated driver helps with precise placement.




These needles are offered in three different sizes (**15mm, 25 mm, 45 mm**) for various uses that can be used in a variety of patient sizes and application sites.



**RT4-2023 –
SUPPLY AND DELIVERY OF EMERGENCY AND RESCUE
EQUIPMENT TO THE STATE FOR THE PERIOD
01 SEPTEMBER 2023 – 31 AUGUST 2028**



ITEM NUMBER	SPECIFICATION	PICTURE
PERSONAL PROTECTION EQUIPMENT		
RT4-01-001	<p><u>SWIFT WATER RESCUE HELMET</u> Helmet for water sport standard: CE En 1385 size: 54/61 cm. Weight: 385 gr. Lifetime: 10 years Technical data sheet :</p> <ul style="list-style-type: none"> • Polystyrene inner shell (eps) • ABS external shell material • Quick release buckle 	 <p style="text-align: center;">CAT NO: RW02</p>
RT4-01-002	<p><u>MOUNTAIN RESCUE HELMET</u> Helmet for mountaineering standard: EN 12492:2012 Size: 54/62 cm. Weight: 350 gr. Lifetime: 10 years Technical data sheet :</p> <ul style="list-style-type: none"> • Polystyrene inner shell (eps) • ABS external shell material • Up/down adjustable size sys Temp for a perfect fitting • Quick release buckle • Chinstrap resistance >50dan • Headlamp clips • Removable chinstrap • Removable and washable inner paddings 	 <p style="text-align: center;">CAT NO: ACCH08.2</p>



ITEM NUMBER	SPECIFICATION	PICTURE
RESCUE EQUIPMENT		
RT4-02-068	<p><u>RESCUE PULLEY (ALUMINIUM) 2.5 -INCH DOUBLE</u></p> <ul style="list-style-type: none"> • Durable and lightweight double pulley for industrial use and hauling systems • For use with static and dynamic textile ropes • Sealed ball bearing and large diameter of the sheave provide great efficiency 94% • Designed and optimized for use with a Prusik friction knot or any clamps • Top attachment point accepts up to three carabiners • Formed side plates provides protection of rope against contact with a surrounding • Fixed-by-forging sheave and anti dust cover of ball bearing provides great lifespan • Side plate holes enable cleaning of the ball-bearing surrounding • Bottom attachment point 12 Kn • Unique production number marking 	 <p>CAT NO: PP08</p>
RT4-02-069	<p><u>RESCUE PULLEY (ALUMINIUM) 2.5 -INCH SINGLE</u></p> <ul style="list-style-type: none"> • Durable and lightweight single pulley for industrial use and hauling systems • For use with static and dynamic textile ropes • Sealed ball bearing and large diameter of the sheave provide great efficiency 94% • Designed and optimized for use with a Prusik friction knot or any clamps • Top attachment point accepts up to three carabiners • Formed side plates provides protection of rope against contact with a surrounding • Fixed-by-forging sheave and anti dust cover of ball bearing provides great lifespan • Side plate holes enable cleaning of the ball-bearing surrounding • Unique production number marking 	 <p>CAT NO: PP07</p>



ITEM NUMBER	SPECIFICATION	PICTURE
RT4-02-070	<p><u>RESCUE PULLEY (ALUMINIUM) 3.5 -INCH SINGLE</u></p> <ul style="list-style-type: none"> • Durable and lightweight single pulley for industrial use and hauling systems • For use with static and dynamic textile ropes • Sealed ball bearing and large diameter of the sheave provide great efficiency 94% • Designed and optimized for use with a Prusik friction knot or any clamps • Top attachment point accepts up to three carabiners • Formed side plates provides protection of rope against contact with a surrounding • Fixed-by-forging sheave and anti dust cover of ball bearing provides great lifespan • Side plate holes enable cleaning of the ball-bearing surrounding • Unique production number marking 	 <p>CAT NO: PP07</p>
RT4-02-097	<p><u>STEEL CABLE SLING – SMALL (3FT)</u></p> <ul style="list-style-type: none"> • 8mm steel cable with PVC coating (total diameter is 8mm) • Workload limit (WLL) – 600kg (safety factor of 6:1) • Soft edge Aluminium ferrules • SANS7531 SWR slings for general purposes • ID tag with Individual Serial number and batch number • Supplied with Batch certificates. 	 <p>CAT NO: PRA02</p>
RT4-02-098	<p><u>STEEL CABLE SLING – MEDIUM (4FT)</u></p> <ul style="list-style-type: none"> • 8mm steel cable with PVC coating (total diameter is 8mm) • Workload limit (WLL) – 600kg (safety factor of 6:1) • Soft edge Aluminium ferrules • SANS7531 SWR slings for general purposes • ID tag with Individual Serial number and batch number • Supplied with Batch certificates. 	 <p>CAT NO: PRA01</p>
RT4-02-099	<p><u>STEEL CABLE SLING – LARGE (5FT)</u></p> <ul style="list-style-type: none"> • 8mm steel cable with PVC coating (total diameter is 8mm) • Workload limit (WLL) – 600kg (safety factor of 6:1) • Soft edge Aluminium ferrules • SANS7531 SWR slings for general purposes • ID tag with Individual Serial number and batch number • Supplied with Batch certificates. 	 <p>CAT NO: PRA03</p>




ITEM NUMBER	SPECIFICATION	PICTURE
RT4-02-100	<p><u>STEEL CABLE SLING – X-LARGE (6FT)</u></p> <ul style="list-style-type: none"> • 8mm steel cable with PVC coating (total diameter is 8mm) • Workload limit (WLL) – 600kg (safety factor of 6:1) • Soft edge Aluminium ferrules • SANS7531 SWR slings for general purposes • ID tag with Individual Serial number and batch number • Supplied with Batch certificates. 	 <p>CAT NO: PRA04</p>
RT4-02-129	<p><u>SPECIALIZED STRETCHER SYSTEM – CONFINED SPACE STRETCHER SYSTEM FOR RESCUE PURPOSES</u></p> <p>The RL1000 is the most basic of the rollup stretchers from Kohlbrat & Bunz. Lightweight, packable and versatile the rollup stretchers are the go-to stretcher for mountain and urban rescue.</p>	 <p>CAT NO: HRSTRL1000</p>
RT4-02-141	<p><u>RESCUE SWIVEL</u></p> <p>Material: Weight: 90g Strength: 22k N Working load: 4kN Color: yellow / antracid aluminum alloy Intended for working at heights and climbing Unique number marking Ball-bearing: working-load 4 k N – Minimum breaking strength: 22 kN</p>	 <p>CAT NO: PR03</p>
RT4-02-142	<p><u>MAILLON</u></p> <p>Maillon Rapide quick links used as a connector on P.P.E. (personal protective equipment) against falls from heights, mountaineering and rock climbing, are distinguished by a specific marking. Breaking load – 8mm: 2800kg / 27.5kn</p>	 <p>CAT NO: PCM01</p>

ITEM NUMBER	SPECIFICATION	PICTURE
RT4-02-171	<p><u>ADJUSTABLE LANYARD FITTED WITH A DIRECTIONAL CONNECTOR WITH WITH VERY LARGE GATE OPENING AND AUTOMATIC LOCKING – LENGTH 3M</u></p> <p>The work positioning lanyard with great handling makes your work safe and efficient.</p> <ul style="list-style-type: none"> • 11 mm rope sewn-in unopenable version of LOCKER • rope clamp - rope covered with tubular webbing • parking cam for blocking in required position • available with stitched-in connector in accordance with EN 362 • Weight: 410 g 560 g (with connector K350) • Strength: 15 kN • Length: 2 m, 4 m • Diameter of rope: 11 mm • Materials: Rock clamp: dual • Rope: polyamide • Colour: Rock clamp: black • Rope: white • Rope protector: black 	 <p>CAT NO: SSA01/4 SCAFF</p>
RT4-02-217	<p><u>RESCUE CARABINERS STAINLESS STEEL AUTO LOCK</u></p> <p>Versatile HMS locking biner designed for belaying, rappelling and anchor building.</p> <ul style="list-style-type: none"> • hot-forged I-beam construction makes the carabiner extremely strong, stiff and light; • large rope bearing area for better rope glide and reduced wear on the carabiner; • pear-shaped basket facilitates belaying with the Munter hitch or other belay devices; • angled snag free Keylock nose for smooth clipping and unclipping manoeuvres; • opening and basket area are large enough to easily accommodate the clove hitch; • twist lock gate prevents from accidental gate opening ; • available also with screw and triple lock gate, all of them easy to handle; • each carabiner is individually tested for the strength 10 kN; • unique production number for better traceability and inspection; • certified for both sport and industrial use. 	 <p>CAT NO: PCLS042</p>

ITEM NUMBER	SPECIFICATION	PICTURE
RT4-02-218	<p><u>RESCUE CARABINERS STAINLESS STEEL SCREW GATE</u></p> <p>Versatile HMS locking biner designed for belaying, rappelling and anchor building</p> <ul style="list-style-type: none"> • hot-forged I-beam construction makes the carabiner extremely strong, stiff and light • large rope bearing area for better rope glide and reduced wear on the carabiner • pear-shaped basket facilitates belaying with the Munter hitch or other belay devices • angled snag free Keylock nose for smooth clipping and unclipping manoeuvres • opening and basket area are large enough to easily accommodate the clove hitch • screw locking sleeve for great handling and effective locking • for better gear management the screw locks are in three different colours • available also with twist and triple lock gate, all of them easy to handle • each carabiner is individually tested for the strength 10 Kn • unique production number for better traceability and inspection • certified for both sport and industrial use 	 <p>CAT NO: PCLS04</p>
RT4-02-220	<p><u>MPD</u></p> <p>The RD2 is a robust self-braking descender designed for experienced users in rescue and rope access.</p> <p>Tested to EN12841:C it has a WLL of 140kg and for rescue purposes a WLL of 200kg. It incorporates a high-efficiency (one-way) pulley wheel with an 88% efficiency on a 1:1 haul. The integrated becket allows the user to set up mechanical advantage systems directly off the device.</p> <ul style="list-style-type: none"> • Weight : 790g • WLL : 140kg or 200kg under rescue conditions • Suitable for 11mm EN1891 rope • Compliance : EN12841:C 	 <p>CAT NO: RD2</p>

ITEM NUMBER	SPECIFICATION	PICTURE
RT4-02-221	<p><u>RESCUE FIGURE OF 8 (ALUMINIUM)</u> A light Alloy special figure eight for the rescue situations and fast rappelling.</p> <ul style="list-style-type: none"> • special shape for multiple rope locking • allows rappelling with heavy loads • high strength • wide range of rope diameters <p>Color: black matt Weight: 154 g • 5.43 oz Material: light-alloy Strength: 40 kN</p>	 <p>CAT NO: PD2</p>
RT4-02-268	<p><u>TWO (2) PERSON TENT</u></p> <ul style="list-style-type: none"> • Sleeps: 2 • Packed weight: 2.35Kg (5lb 3oz) • Packed size: 32cm x 21cm • Season: 3 season • Number of porches: 2 • Number of doors: 2 • Flysheet: 68d 210T R/S PU 4000mm FR WR C6 • Inner: 68d 210T R/S BR / FR • Groundsheet: 68d 190t PU 5000mm FR WR C6 • Poles: 8.5mm Superflex Alloy • Pegs: 10 x Alloy 'V' Angle • Guylines: 4 x Yellow Reflective • Pitch type: Fly and inner pitch together 	 <p>CAT NO: TDW36.01</p>

ITEM NUMBER	SPECIFICATION	PICTURE
RT4-02-270	<p><u>SELF-INFLATING HIKING MATTRESS</u> Light and compact self inflating mattress providing desirable comfort and support as well as favourable packed size. Ranges from 2cm/0.75' to 3.8cm/1.5' thickness. Various shapes suit different needs. For Rectangular model, Diamond die cut for weight saving. For Mummy model, ergonomic designed shape fits your body as well as your sleeping bag. Dot die cut, tapered top and bottom help to shave extra weight. Anti-fungal treated foam.</p> <ul style="list-style-type: none"> • Inflates to 3.8cm thickness • Screw valve • Mummy shaped 75D nylon • Diamond Ripstop fabric Diamond die cut foam reduces weight while retaining insulation & support • Anti-fungal treated foam • Repair kit included • Stuff bag included • Velcro tie strap included • Estimated R-Value: ~ 3.8 (1.5/3.8cm) Size: 183cm Long x 51cm Wide x 3.8cm Thick • Pack size: 29 x Ø18cm • Weight: 795g • Colour Red 	 <p>CAT NO: MAT-LMT015</p>
RT4-02-274	<p><u>ROPE BAG – 100M</u></p> <ul style="list-style-type: none"> • A large haul bag with shoulder straps and haul loop. • A simple closing mechanism, drainage hole and brightly coloured in either orange or yellow. 	 <p>CAT NO: RB003.2</p>

ITEM NUMBER	SPECIFICATION	PICTURE
RT4-02-286	<p><u>ENERGY ABSORBING LANYARD</u> Lightweight, practical shape does not restrict users at the work. Designed for users with a maximum weight of 140 kg (including equipment). Optional length of lanyards for diverse work. Newly tested for horizontal use and a resulting simulated fall over edge</p>	 <p>CAT NO: SLSA02.1</p>
RT4-02-287	<p><u>25MM LIFE SAFE WEBBING</u> The TUBULAR TAPE is a 26mm wide American type strap. This system allows for soft, flexible straps with a good grip.</p> <ul style="list-style-type: none"> Standards : CE EN565 Reference : BST26 Material : Polyamide Weight per meter : 38g Breaking load: 1500daN Colours: Red, Blue, Black 	 <p>CAT NO: SW26/2</p>
RT4-02-289	<p><u>RESCUE HIKING BACKPACK</u> The Vango Sherpa 70:80 rucksack is one of the most popular hiking rucksacks. Updated for 2022, Sherpa joins Vango's Earth Collection. The Sherpa is extremely practical and durable. Providing optimum fit through its Self-Adjust Back System and a steady air flow through the shoulder straps due to Air Mesh, it's the perfect companion for those doing multi-day treks and activities. The Sherpa is part of the recommended kit by The Duke of Edinburgh's Award. The Vango Hydrant Pack is the perfect accessory to your hydration compatible bag</p>	 <p>CAT NO: BPKSHERP.01</p>

Kamo, Bongeka, Cami	Internal Sales	sales@statmedical.co.za	011 708 6363 066 186 4334
Dagan Warns	National Capex Manager	dagan@statmedical.co.za	083 600 5990
Sihle Biyela	KZN Sales Consultant	kzn@statmedical.co.za	079 832 4085
Leon Swanepoel	General Manager	leon@statmedical.co.za	083 230 7322
Stanley Trott	Clinical Product Specialist EMS	ems@statmedical.co.za	083 601 6298
Doreen Timmerman	Customer Service Manager	doreen@statmedical.co.za	083 600 3737
Justin Bannister	Sales Consultant	sales2@statmedical.co.za	083 605 5851
Mokgadi Sebata	Sales Technician	sale1@statmedical.co.za	083 679 3622
Joseph Soga	Sales Consultant Western Cape	wcsales@statmedical.co.za	083 600 7521
Wiseman Ngcebesha	Warehouse Manager	wiseman@statmedical.co.za	083 600 3741
Cathy Crawford - Nutt	National Disposables Manager	cathy@statmedical.co.za	083 391 1155



227 Main Avenue Ferndale,
Randburg, 2194
T: +27 11 708 6363