



OXYLITRE PSP001 "petite elite" Portable Suction Pump High Suction STANDARD Re-Usable Jar



Portable Suction Pumps

Oxylitre's range of Portable Suction Pumps (PSP's) deliver a high quality flexible portable suction solution for Healthcare and Emergency Services. Our contemporary design Portable Suction is user friendly, comprehensively field tested to the latest standards. Made from high impact plastic, the units streamline compact design is both durable and lightweight. The Portable Suction Pumps (PSP's) are intuitive to use incorporating a luminous LED illuminated touch control panel and an easy access battery compartment. These units are reliable and extremely functional our range of Portable Suction Units offer high performance with operational reliability.

Oxylitre's range of Portable Suction Units are dependable, with minimum maintenance required. Our Range of Portable Suction Units come with a long service life and a 7-year warranty.



Specification

High Vacuum High Flow
Single 1000ml Standard Autoclavable Jar Assembly
Robust High Impact Resistant Moulded Enclosure
Encapsulated Bacterial/Hydrophobic Filter Disk
Luminous Touch Panel Control
Instant Maximum Vacuum Override Button
Battery Level/ Charging Indicator
80 to 500mmHg.
Flowrate >25Lpm
110volt – 240volt AC
12volt – 28volt DC
12 volt DC Rechargeable Battery, 4.3mAh
Run time >1hr 20min
Weight 4.9Kg
Approximate Size 310mm x 305mm x 170mm
Peak Sound Output at MAX Vacuum from 1m: 65 decibels

STAT-TIAKENI MEDICAL PTY LTD

Reg. No 2007/010669/07 | VAT No. 4090241185

227 Main Avenue, Ferndale, 2194

P + 27 (0) 11 708 6363 | F +27 (0)11 708 6364 | E sales@statmedical.co.za

Directors: SJ McDonnell, BW Poynton, P Fourie, G Leonard

